

NAME

peakdeck.pl v1.1 - a kernel density estimator based peak caller for DNaseI-seq data.

**SYNOPSIS****Numerical sorting:**

```
peakdeck.pl -NS mappedData.sam > sortedData.sam
```

Filtering:

```
peakdeck.pl -F sortedData.sam -g hg19.chrom.sizes -q 15 -PCR ON > filteredData.sam
```

Random read selection:

```
peakdeck.pl -R filteredData.sam -nr 20000000 > filteredData.20m.sam
```

Density analysis:

```
peakdeck.pl -D filteredData.20m.sam -g hg19.chrom.sizes > densityTrack.wig
```

Peak calling:

```
peakdeck.pl -P filteredData.20m.sam -g hg19.chrom.sizes > peakList.bed
```

Size ordering of peaks:

```
peakdeck.pl -T peakList.bed -n 50000 > orderedPeakList.bed
```

ARGUMENTS**Numerical sorting (-NS mappedReads.sam)**

Sorts sam format reads by base start position, irrespective of chromosome. Equivalent to the Linux/Unix/osx command: [sort -n --key=4,5 filename.sam > sortedFilename.sam]. For fast results, memory corresponding to ~2.5 times file size should be available.

Sam filtering (-F mappedReads.sam)

Sorts sam files by chromosome, in the order that chromosomes appear in the chromosome size file. The chromosome size file is mandatory. The chromosome size file is a plain text, tab separated file in the format:

chr1	249250621
chr2	243199373
chr3	198022430
.....
chrN	size(bp)

Mandatory settings

-g /path/to/chromosomeSizeFile.txt

Specifies the path to the text file containing tab-separated list of chromosome names and sizes.

Optional settings

-q integerValue

Specifies a mapq cutoff score for filtering. Reads with a mapq score less than the supplied value will be removed from the resulting filtered file. By default, -q is set to zero, so no filtering for mapq scores will occur.

-u integerValue

Specifies a UQ base mismatch score for filtering. Reads with mismatch scores greater than this value will be removed from the filtered dataset. By default, -u is set to 10000, so that no filtering by uq score will occur.

-i samHeaderFile.sam

Specifies a file containing a sam header, which if set, will be included at the beginning of the newly filtered file.

-PCR ON|OFF

Allows PCR duplicate reads to be removed from sam file. Reads are considered PCR duplicates if adjacent reads have identical chromosome, start position, mapq score, and

sequence. By default -PCR is set to OFF, so no filtering of PCR duplicates will occur. To detect PCR duplicates, chromosomes must be in numerical order (see Numerical sorting above).

Random read selection (-R mappedReads.sam)

Randomly selects a target number of reads from a specified sam file. Selected reads are printed to STDOUT by default.

Mandatory settings

- nr integer
Specifies the target number of reads to be randomly selected from the given sam file. The number of reads must be a positive whole number.

Density analyzer (-D mappedOrderedReads.sam)

Creates a smoothed, unitless read density track in wig format, representing the distribution of reads in the given sam file. Sam files must be grouped by chromosome, and ordered by read start position (see Numerical sorting and Sam filtering above). The order of chromosomes in the density track is determined by the order in which they appear in the mandatory chromosome size file (see Sam filtering for chromosome size file format). By default, the results are printed to STDOUT.

Mandatory settings

- g /path/to/chromosomeSizeFile.txt
Specifies the path to the text file containing tab-separated list of chromosome names and sizes.

Optional settings

- n positiveInteger
Specifies the one-tailed size of the smoothing bin. By default, -t is set to 150, giving a bin size of 300 bp. This value determines both the size of sampling bin, and the width of the Gaussian probability density function used to calculate read densities, and must be a positive whole number.
- STEP positiveInteger
Specifies the size of steps by which the probability density function and sampling bin move along the genome. By default, -STEP is set to 100. Smaller step sizes proportionately increase the number of calculations carried out, and therefore the time taken for the analysis. -STEP must be a positive whole number.
- d positiveInteger
Specifies the standard deviation of the probability density function. This value determines how broadly the read density scores are spread over each sampling bin, and therefore determines the degree of smoothing that occurs. By default -d is set to 50, and must be a positive whole number.
- t positiveInteger
Specifies a low threshold, below which read density scores won't be included in probability density function calculations. By default, -t is set to the number of reads expected to occur in the set bin size if the number of reads in the dataset were randomly distributed. All reads present in the data set will be included in the analysis if -t is set to 0. -t must be a non-negative whole number.
- m positiveInteger
Specifies a high threshold, above which read density scores won't be included in probability density function calculations. By default, -m is set to 100000000, ensuring that no reads will be excluded from analysis in default settings.
- o integer
Specifies a track offset. All positions in the resulting wig file will be offset by this value. For DNaseI-seq data, the read start sites are considered DNaseI cutting sites, and so by default, -o is set to 0. If the centre of the DNA fragment is considered the point of interest (for example, in ChIP-seq), setting -o to half the average fragment size may give a more precise depiction of signal localisation.

Peak calling (-P mappedOrderedReads.sam)

Identifies peaks in the provided sam file, and provides output in bed format to STDOUT. Sam files must be grouped by chromosomes, and ordered by read position (see Numerical sorting and Sam filtering above). The order of chromosomes in the peak file is determined by the order in which they appear in the mandatory chromosome size file (see Sam filtering for chromosome size file format).

Mandatory settings

-g /path/to/chromosomeSizeFile.txt
Specifies the path to the text file containing tab-separated list of chromosome names and sizes.

Optional settings

-bin positiveInteger
Specifies the size of the central sampling bin. By default, -bin is set to 300, which represents the expected average feature size. -bin must be set to a positive whole number

-back positiveInteger
Specifies the size of the background sampling bin. By default, -back is set to 3000, ten times the size of the central sampling bin. -back must be set to a positive whole number and must be larger than the size of the central bin.

-STEP positiveInteger
Specifies the size of steps by which the sampling bin moves along the genome. By default, -STEP is set to 100. Smaller step sizes proportionately increase the number of calculations carried out, and therefore the time taken for the analysis. -STEP must be a positive whole number.

-FLAT positiveInteger
Specifies a flat threshold for peak calling in reads per bin. When -FLAT is set, the threshold calculated by PeakKDEck for peak calling is overridden, and the value given by -FLAT is used in its place. FLAT must be a positive number.

-b /path/to/blueprintFile.bed
This option provides the path to a bed file which contains a list of contiguous genomic loci indicating the sites of known open chromatin sites, tagged with the number of cell types with open chromatin at that site. The format is as follows:

chr1	1	10099	C#1	0
chr1	10100	10330	#1	37
chr1	10331	10344	C#2	0
chr1	10345	10590	#2	4
chr1	10591	16099	C#3	0

where the columns respectively indicate the chromosome name, site start position, site end position, element name, and number of cell types with open chromatin at that site. When this file is provided, PeakKDEck calculates signal-to-noise ratio, and calculates the background probability distribution from sites selected from loci with no known open chromatin.

-npBack positiveInteger
Sets the number of sites to randomly select to calculate the background probability distribution. By default this is set to 50000 sites. -npBack must be a positive whole number.

-sig probabilityValue
Specifies the positive limit of the probability distribution for selecting the corrected read density for peak threshold. By default, -sig is set to 0.001. -sig must be a positive number between 0 and 1.

-PVAL ON|OFF
Peaks are scored with the maximum corrected read density recorded during that peak by default. Setting -PVAL to ON converts this corrected read density to a probability value from the background probability distribution used to calculate the threshold. This value represents the probability that a corrected read density of that magnitude belongs to the background probability distribution.

Top peak selection (-T peaks.bed)

Sorts peak bed files in descending order by corrected read density score. By default, the sorted peaks are printed in bed format to STDOUT. The target file must be in bed format.

Mandatory settings

Optional settings

-n positiveInteger

This specifies the number of peaks to include in the resulting bed file, from the highest scoring peak downwards. By default, -n is set to ALL, and all the peaks are printed to the output file. -n must be either 'ALL' or a positive whole number.

DESCRIPTION

PeakKDEck is a utility written in Perl, mainly intended for use in the identification of peaks in mapped DNase-seq data. It also includes a set of utilities for processing and manipulation of this data from the mapping stage forwards. It works on data in sam format.

PeakKDEck selects a threshold read density for peak calling by constructing a probability distribution of background read density scores using kernel density estimation. It selects a threshold by selecting a read density that is 'significantly' outside this background probability distribution. All measurements of read density are corrected for local background variation in signal intensity.

PeakKDEck is also available as a standalone GUI application for all major platforms. The GUI wrapper is written in Perl/Tk, and is available at www.ccmp.ox.ac.uk/PeakKDEck.

FAQs

What are the system requirements?

The command line and GUI PeakKDEck applications have been tested on OSX (Mountain Lion), Ubuntu 12.04 LTS, Windows XP and Windows 7. The system requirements are largely dependent on the size of data files being used. We recommend at least 4GB memory for basic use with small data files. For the numerical sorting of sam files, $\sim(\text{file size} * 2.5)$ free memory is required for efficient sorting. For the command line applications, we recommend Perl v5.12 or later. On Windows, PeakKDEck was tested with Strawberry Perl.

How do I install PeakKDEck?

PeakKDEck GUI: On Windows, PeakKDEck should run without the need to install Perl, or other additional software's. On Linux and OSX platforms, the X Window System (X11 or XQuartz) must be installed.

PeakKDEck command line: to use the PeakKDEck command line application, Perl must be installed on your computer. We recommend Perl v5.12 or later. On Linux and OSX platforms, no other software is required to run the command line application. On the Windows platform, a pseudorandom number generating module (Math::Random::MT) is required, and is available through CPAN for Strawberry Perl users, and PPM (Math-Random-MT) for ActiveState users.

Which short read file formats does PeakKDEck work with?

At present, PeakKDEck only works with files in the SAM format (see samtools.sourceforge.net/SAMv1.pdf for details).

Which application should PeakKDEck be opened with?

On the OSX platform, after launching the PeakKDEck GUI application (having installed XQuartz), you may be prompted to choose an application with which to open PeakKDEck. PeakKDEck should be opened with the Terminal application (located at /Applications/Utilities/Terminal.app).

Why is the GUI freezing?

As yet, the PeakKDEck GUI is not a multi-threaded program. As such during data processing, the GUI may appear frozen or unresponsive. Particularly on the Windows platform. For now, this is expected behaviour.

The GUI will refresh when new status updates are available, and will return to full responsiveness when data processing has finished.

EXAMPLES

see SYNOPSIS above.

CAVEATS

AUTHOR

Michael McCarthy
michaeltmccarthy@gmail.com

ACKNOWLEDGEMENTS

SEE ALSO